



< Example of DNA extraction >



< Kaneka Easy DNA Extraction Kit version 2\* >



- Simple protocol
- 10 minute extraction

\*This product is for research use only.


<Performance comparison with a spin column DNA purification kit>

Extract DNA from HEK293 cells by using Kaneka Easy DNA Extraction Kit version 2\* following PCR.

Target sequence: endogenous gene (CYBB E2, E5, E10)



Electrophoresis for Confirming DNA amplification



Kaneka DNA Easy Extraction Kit version 2 successfully extracted DNA from HEK293 cells..